


Ucha05 Ucha10 – Simple chargers for 10-150 mA cells

These low cost single cell charger is powered by any USB port. A connector can be added to receive the current from a 5 to 6V DC supply.

Bahoma magnets charges 5mm and 10mm Plantraco Lipos. A connector or direct wires can be soldered if you do not use Bahoma.

Ucha05 uses 5mm pitch Bahoma and is planned for 20 to 50mA Lipo. Charging current is fixed, 35mA and can be modified by unsoldering a resistor.

Ucha10 uses 10mm pitch Bahoma and is planned for 80 to 150mA Lipo. Charging current is fixed, 100mA. USB is limited to 100 mA, it is not recommended to change the resistor. The charger circuit is limited to 150 mA if an external supply is used.



Price 20 CHF ca 20 USD

Quick Lipo tester

Drop your Lipo Both Leds On : voltage > 3.8V

5mm Bahoma for 20-50 mAh cells

10mm Bahoma for 70-200 mAh cells

Press right switch and discharge with 22 Ohm (160 mA)
 Press left switch and charge with 10 Ohm (330 mA)



Press right switch and discharge with 12 Ω (300 mA)
 Press left switch and discharge with 6 Ω (500 mA)

Price 12 CHF ca 12 USD

Other products by Didel for less than 150 mAh Lipos

Bicha



Price 38 CHF ca 38 USD

Bapi-Bami pieces



Price 1.30